


ANTERION® 生物测量

检查质量评估指南



Anterior Axial Curvature (3 mm Ring)

R _{mean} ¹⁻⁷	8.02 – 7.57 mm
SimK _{mean} ¹⁻⁷	42.1 – 44.6 D

0.17 mm 的测量误差将在正常眼中导致 1 D 术后误差。

双眼之间 SimK 差值 > 1 D 并不常见。²⁶

1 D SimK 测量误差在正常眼中导致 1 D 术后误差。

Posterior Axial Curvature (3 mm Ring)

K _{mean} ¹⁻¹⁰	-5.6 – -6.6 D
P/A ratio ²⁰	0.82 – 0.86

P/A 比率 可通过角膜屈光手术改变。该比率将在近视治疗后降低,而在远视治疗后增高。

Total Corneal Wavefront (6 mm)

Sph. aberration ²¹	0.25 – 0.27 μ m
RMS HOA ²¹	0.45 – 0.48 μ m

Pachymetry

CCT (vertex) ^{1-4, 6-15}	474 – 608 μ m
-----------------------------------	-------------------

Anterior Segment

AQD (ACD) ^{4,15-19}	2.11 – 3.91 mm
从角膜后表面到晶状体前表面	
Lens thickness ^{4, 12-15, 22}	3.43 – 4.77 mm

晶状体厚度随年龄增长而增加。²⁷

Axial Length

Length ^{4,10-15,22}	21.5 – 26.4 mm
------------------------------	----------------


眼轴长度短出现于轴性远视中,但也可能由例如视网膜色素上皮 (RPE) 隆起引起。轴性近视以眼轴长度长为特征。

White-to-white


WTW ^{2-4, 12-14, 22-25}	11.0 – 12.8 mm
----------------------------------	----------------

双眼之间眼轴长度差值 > 0.5 mm 并不常见。²⁶ 检查是否存在屈光参差。

1 mm 测量误差将在正常眼中导致 2.7 D 术后误差。



Axial length: **23.52 mm** $\pm$ 0.002 mm
Eye status: No surgery, Phakic, Vitreous only



参考文献

1. KIM KY et al. Anterior segment characteristics in normal and keratoconus eyes evaluated with a new type of swept-source optical coherence tomography.PLoS One 2022; 17(9):e0274071.DOI:10.1371/journal.pone.0274071.
2. HERBER R et al. Agreement and repeatability of corneal tomography in healthy eyes using a new swept-source OCT a rotating Scheimpflug camera and a dual Scheimpflug-Placido system.J Cataract Refract Surg 2022; 48(2):190-198.
3. KIM K et al. Diagnostic Validation of the Screening Corneal Objective Risk of Ectasia Analyzer Evaluated by Swept Source Optical Coherence Tomography for Keratoconus in an Asian Population.Bioengineering (Basel) 2023; 10(11):1335.
4. The Heidelberg Engineering ANTERION Anterior Segment Cornea and IOL Precision and Agreement Study [data on file].
5. TAÑÁ-RIVERO P et al. Repeatability of whole-cornea measurements using a new swept-source optical coherence tomographer.Eur J Ophthalmol 2021; 31(4):1709-1719.
6. ESCOLANO SERRANO J et al. Intraobserver Repeatability of Tomographic Pachymetric and Anatomical Measurements in Healthy Eyes Using a New Swept-Source Optical Coherence Topographer.Cornea 2022; 41(5):598-603.
7. HERBER R et al. Comparison of corneal tomography using a novel swept-source optical coherence tomographer and rotating Scheimpflug system in normal and keratoconus eyes: repeatability and agreement analysis.Eye Vis (Lond) 2022; 9(1):19.
8. PÉREZ-BARTOLOMÉ F et al. Anterior-Segment Swept-Source Ocular Coherence Tomography and Scheimpflug Imaging Agreement for Keratometry and Pupil Measurements in Healthy Eyes.J Clin Med 2021; 10(24):5789.
9. SAAD A et al. Discrimination between keratoconus forme fruste keratoconus and normal eyes using a novel OCT-based tomographer.J Cataract Refract Surg 2023; 49(11):1092-1097.
10. SCHIANO-LOMORIELLO D et al. Repeatability of automated measurements by a new anterior segment optical coherence tomographer and biometer and agreement with standard devices.Sci Rep 2021; 11(1):983.
11. CHENG SM et al. Repeatability of a new swept-source optical coherence tomographer and agreement with other three optical biometers.Graefes Arch Clin Exp Ophthalmol 2022; 260(7):2271-2281.
12. TAÑÁ-RIVERO P et al. Agreement between 2 swept-source OCT biometers and a Scheimpflug partial coherence interferometer.J Cataract Refract Surg 2021; 47(4):488-495.
13. TAÑÁ-RIVERO P et al. Lens-vault analysis and its correlation with other biometric parameters using swept-source OCT.J Optom 2022; 15(1):88-99.
14. TANA-SANZ P et al. Agreement of predicted intraocular lens power using swept-source optical coherence tomography and partial coherence interferometry.Expert Rev Med Devices 2021; 18(12):1219-1234.
15. RUIZ-MESA R et al. Ocular biometric repeatability using a new high-resolution swept-source optical coherence tomographer.Expert Rev Med Devices 2020; 17(6):591-597.
16. CHENG SM et al. Repeatability and Agreement of Two Swept-Source Optical Coherence Tomographers for Anterior Segment Parameter Measurements.J Glaucoma 2022; 31(7):602-608.
17. XIE X et al. Age- and refraction-related changes in anterior segment anatomical structures measured by swept-source anterior segment OCT.PLoS One 2020; 15(10):e0240110.DOI:10.1371/journal.pone.0240110.
18. PARDESHI AA et al. Intradevice Repeatability and Interdevice Agreement of Ocular Biometric Measurements:A Comparison of Two Swept-Source Anterior Segment OCT Devices.Transl Vis Sci Technol 2020; 9(9):14.
19. CHAN PP et al. Anterior chamber angle imaging with swept-source optical coherence tomography:Comparison between CASIAII and ANTERION.Sci Rep 2020; 10(1):18771.
20. SAVINI G Posterior Corneal Measurements by 4 Devices.IOL Power Club; November 2023; Palm Springs (not published).
21. PÉREZ-BARTOLOMÉ F et al. Agreement between anterior segment swept source-OCT and Scheimpflug imaging corneal aberration measurements in healthy eyes.Eur J Ophthalmol 2022; 32(6):3363-3371.
22. YUN JS et al. Evaluation of angle-to-angle and spur-to-spur using swept source optical coherence tomography in different refractive error.PLoS One 2022; 17(11):e0277703.
23. MONTÉS-MICÓ R et al. Angle-to-angle and spur-to-spur distance analysis with high-resolution optical coherence tomography.Eye Vis (Lond) 2020; 7:42.
24. MONTÉS-MICÓ R et al. Assessment of anterior segment measurements using a high-resolution imaging device.Expert Rev Med Devices 2020; 17(9):969-979.
25. TAÑÁ-RIVERO P et al. Agreement of white-to-white measurements with swept-source OCT Scheimpflug and color LED devices.Int Ophthalmol 2021; 41(1):57-65.
26. YI L & FANG JB Interocular symmetry analysis of bilateral eyes. Journal of Medical Engineering & Technology 2014; 38(4):179-187.
27. HOFFER, KJ Axial dimension of the human cataractous lens.Archives of Ophthalmology 1993; 111(7):914-918.